

# iCARE Ovulation (LH) Test Kit

JAL [jalinnovation.com/icare-menopause-fsh-cassette-test-kit/](http://jalinnovation.com/icare-menopause-fsh-cassette-test-kit/)

25/02/2013

**Category Name :** Others Rapid Screen Test Kit

The iCARE One Step Ovulation Urine Test is a qualitative test used to predict when there is a LH surge, and in turn, when you are likely to ovulate.

## HOW DOES IT WORK?

The LH (Luteinizing Hormone) which is in the urine of normal women will increase dramatically in the middle of the menstrual cycle. The LH increase triggers ovulation which is when the egg is released periodically from Normal fertile women. WHO experts state that LH testing is a reliable way to detect ovulation. Ovulation will occur after 24-48 hours following a positive test.

When specimen is dropped into cassette, capillary action carries the specimen to migrate along the membrane. When LH in the specimen reaches the test zone region of the membrane, it will form a colored line. Absence of this colored line suggests a negative result.

To serve as a procedure control, a colored line will appear at the control zone region, if the test has been performed properly.

The test device detects ovulation with a high degree of certainty. It is a valuable tool in helping to achieve pregnancy since it determines the period that the egg and the sperm cells will meet in the best conditions.

Ovulation may be irregular because of circumstances, emotions and other factors in your life. You cannot presume that ovulation always occurs at the same time after menstruation. Therefore, you should test again in each menstrual cycle.



## Specifications

Category Name	Others Rapid Test Kits
Item Number	WF1315 (Cassette) WF1314 (Midstream)
Sample	Urine Only
Format	Test card
Accuracy	99%
Shelf Life	24 months from date of production

## Specifications

---

Storage 2 to 30 C in sealed pouch under dry conditions

---

Certifications FDA510k CE1023

---

Made by [JAL Innovation \(S\) Pte Ltd](#)