

iCARE ABO & RhD Blood Grouping Test Card

JAL jalinnovation.com/icare-abo-rhd-blood-grouping-test-card/

01/03/2013

Category Name : Others Rapid Screen Test Kit

The assay starts with a sample applied to the sample well and added sample diluents provided. The monoclonal anti-A and anti-B immobilized respectively on the pad in sample region can react with the corresponding antigen on the surface of the red blood cells (RBC). A positive sample produces an immune response and the RBC can be captured on the sample region as a red signal, indicating that the test is positive. A negative sample does not produce an immune response, after adding sample diluents, the RBC can be washed away and the absence of the RBC indicates that the test is negative.



Features

Short Reaction Time

Easy to Read Results

High Sensitivity 100%

No Required Aids and Equipment

Room Temperature Storage 4-30 C

Low Risk of Pollution

Shelf Life 24 Months

Low Cost

Specifications

Category Name	Others Rapid Test Kits
Item Number	IT1064
Sample	Whole blood / Serum / Plasma
Format	Test card
Type of Test	Solid Phase
Shelf Life	24 months from date of production
Storage	2 to 30 C in sealed pouch under dry conditions

Specifications

Made by [JAL Innovation \(S\) Pte Ltd](#)