

OK Match II Blood Glucose Meter

JAL jalinnovation.com/ok-match-ii-blood-glucose-meter-2/

25/02/2013

Features No coding Alternative Site Testing

Accurate result in 6

seconds Light weight Very small blood sample, only 0.7

ul 250 test memories with date and time Large display

Risk indicator for lower than 70 mg/dL or

higher than 180 mg/dL 5 years' warranty PC download



System Specifications

Model Name	OK Match II
Measurement Range	20 – 600 mg/dL (1.1 – 33.3 mmol/L)
Specimen Type	Capillary Whole Blood
Alternative Site Testing	YES (palm, forearm, upper arm, calf, or thigh)
Calibration Against	Plasma-equivalent
Measuring Time	6 seconds
Sample Volume	0.7 ul
Sensitivity in HCT Range	20-60%
Operating Conditions	10 C to 40 C (50 F to 104 F), between 10 to 85% R. H. Altitude ≤ 3,402 m (11,161 ft) above sea level
Storage / Delivery Conditions	10 C to 40 C (50 F to 104 F), between 10 to 85% R. H. (Strip included). 20 – 50 C (4 – 122 F) ≤ 90% R.H. (Strip excluded)
Memory Capacity	250 test results with time and date

System Specifications

Risk Indicator	<ul style="list-style-type: none">– When test within 70-120 mg/dL (3.7 – 6.9 mmol/L)– When test \leq 70 mg/dL (3.9 mmol/L) or $>$ 180 mg/dL (10 mmol/L)
Average Calculation	7, 14, 21, 28 days
Downloadable	YES (RS-232 port)
Battery Type	3V Li-battery (CR 2032) X1
Battery Life	1,000 times of measurements
Auto OFF	3 minutes without action
Meter Dimension	99mm(L) x 62mm(W) x 20mm(H)
Meter Weight	Approximate 50g (with battery)
Expiration for test strip	Storage : 24 months In-Use 180 days

[Previous](#)[Next](#)



Made by [JAL Innovation \(S\) Pte Ltd](#)