

iCARE Gonorrhoea Rapid Screen Test

JAL jalinnovation.com/icare-gonorrhoea-rapid-screen-test/

Category Name : Sexually Transmitted Disease (STD) Rapid Screen Test Kit

The iCARE One Step Gonorrhoea rapid home Screen test (swab) is a rapid home test kit for direct detection of Gonorrhoea.

This product is for in vitro diagnostic use only

Gonorrhoea is a common human sexually transmitted infection caused by the bacterium *Neisseria gonorrhoeae*.

With our iCARE Gonorrhoea Rapid Screen Test kit you can detect Gonorrhoea on your own in 100% privacy. Our test kit is easy to use and the result can be read within 15 minutes



Specifications

| | |
|---------------|--|
| Category Name | STD Rapid Test Kit |
| Catalog No | WF1002 |
| Packing | Hospital version: 50T Field version: 1T OTC Packing |
| Specimen | Whole blood / Serum / Plasma |
| Type of Test | Antigen |
| Format | Hospital version: 1 well cassette, pipette, diluents (4 bottles) Field version: Buffer A, Buffer B, extraction tube, swab, 1 well cassette, pipette |
| Sample | 10ul serum |
| Result | Within 5 minutes |
| Specificity | 98% |
| Sensitivity | 98.8% |

Specifications

Shelf Life 24 months

Storage Store at 4 °C ~ 30 °C in the sealed
pouch up to the expiration date.

Made by [JAL Innovation \(S\) Pte Ltd](#)