



iCare

iCARE RAPID IMMUNOGRAPHIC ASSAY MULTI-INFECTIOUS COMBO TEST (HIV 1&2 . HBSAG . HCV. TP)

Multi-Infectious Disease Test (HIV/HBsAg/HCV/TP)

Rapid, Immunochromatographic Assay for Qualitative Detection of Infectious Disease (HIV/HBsAg/HCV/TP) in human whole blood, serum or plasma

Human Immunodeficiency Virus (HIV)

- Sensitivity: 100%
- Specificity: 100%
- Accurate results in 15 minutes

Herpatitis B (HBsAg)

- Sensitivity: 100%
- Specificity: 99%
- Accurate results in 15 minutes

Herpatitis C (HCV)

- Sensitivity: 100%
- Specificity: >99%
- Accurate results in 15 minutes

Syphilis (TP)

- Sensitivity: 100%
- Specificity: 99%
- Accurate results in 15 minutes

Test Procedures



Apply sample to sample well




Add diluent to diluent well



Read results in 15 minutes

Interpreting Results



Positive

Negative

Invalid

Smart Options For Whatever You Want ...



Blood donation screening



Antenatal care system



Routine examination

